

Date: 2026-08-07
Request ID: 100048806

Microwave				
Analyte(s)	Rh			
SOP#				
Acid Cocktail	6 mL HCl, 2 mL HNO3			
Rack	7			
Vessel	Quartz			
Stir	90%			
Time (min)	Power (W)	T1 (°C)	T2 (°C)	P (bar)
15	1200	180	65	110
15	1500	240	65	110
30	1500	240	65	110
sample(s)	weight (g)	volume (mL)		
200128690_1	0.3081	50		
200128690_2	0.3071	50		

Date: 2026-08-07
Request ID: 100048806

LOI temp: 900 °C

sample	crucible (g)	crucible + sample (g)	after 900 °C	LOI (%)	correction
200128690	17.9925	19.9928	19.6521	17.03	0.8297

Rh ICP analysis

sample	Dilution	conc. [mg/L]	Rh_ppm [dried]	%Rh [dried]
200128690_1	10/3	1.332	868.47	0.09
200128690_2	10/3	1.361	890.26	0.09

rh result: 879.37 ppm

rh result: 0.09 %

rh lot: *Manufacturer: Ricca, 45120966 exp: 2027-11-30*

Note(s):

$$\text{ppm } M+ = [\text{conc.}][\text{volume}][\text{dilution}]/[\text{weight}]$$

$$\%M+ = [\text{conc.}][\text{volume}][\text{dilution}]/([\text{weight}]*10,000)$$

$$\%MO = [\text{oxide factor}]*\%M+$$

A = crucible

B = crucible + sample

C = crucible + sample after drying

$$\%LOI = ([B-A]-[C-A])/(B-A)*100\%$$

$$\text{ppm } M+ \text{ [dried]} = \text{ppm } M+ / ([1-(\%LOI)/100])$$

$$\%M+ \text{ [dried]} = \%M+ / ([1-(\%LOI)/100])$$

$$\%MO \text{ [dried]} = \%MO / ([1-(\%LOI)/100])$$