

Date: 2026-05-27
Request ID: 100048585

Microwave				
Element(s)	Ca			
SOP#				
Acid Cocktail	6 mL HCl, 2 mL HNO3			
Rack	7			
Vessel	PTFE			
Stir	90%			
Time (min)	Power (W)	T1 (°C)	T2 (°C)	P (bar)
15	1200	180	65	110
15	1500	240	65	110
30	1500	240	65	110
sample(s)	weight (g)	volume (mL)		
200128098_1	0.1101	250		
200128098_2	0.1021	250		
200128099_1	0.1044	250		
200128099_2	0.1081	250		

Date: 2026-05-27
Request ID: 100048585

Ca ICP analysis

sample	Dilution	conc. [mg/L]	Ca_ppm	%Ca
200128098_1	100/3	1.950	147593.10	14.76
200128098_2	100/3	1.799	146833.17	14.68

ca result: #####

ca result: 14.72 %

ca oxide factor: 1.399

ca oxide result: 20.60 %

ca lot: *Manufacturer: Assurance, 28-35CAY exp: 2026-10-30*

Note(s):

Date: 2026-05-27
Request ID: 100048585

Ca ICP analysis

sample	Dilution	conc. [mg/L]	Ca_ppm	%Ca
200128099_1	100/3	1.823	145514.05	14.55
200128099_2	100/3	1.857	143154.49	14.32

ca result: #####

ca result: 14.43 %

ca oxide factor: 1.399

ca oxide result: 20.19 %

ca lot: *Manufacturer: Assurance, 28-35CAY exp: 2026-10-30*

Note(s):

A = crucible

B = crucible + sample

C = crucible + sample after drying

$$\%LOI = ([B-A] - [C-A]) / (B-A) * 100\%$$

$$ppm\ M+ = [conc.][volume][dilution] / [weight]$$

$$\%M+ = [conc.][volume][dilution] / ([weight] * 10,000)$$

$$\%MO = [oxide\ factor] * \%M+$$

$$ppm\ M+ \ [dried] = ppm\ M+ / ([1 - (\%LOI) / 100])$$

$$\%M+ \ [dried] = \%M+ / ([1 - (\%LOI) / 100])$$

$$\%MO \ [dried] = \%MO / ([1 - (\%LOI) / 100])$$

$$\%Cr(VI) = (1.733[mL\ FAS][N\ FAS]) / ([weight])$$

$$\%Cr_2O_3 = 1.462 * [\%Cr(VI)]$$

$$\%Cr(III) = total_Cr - Cr(VI)$$