

Date: 2026-06-22

Request ID: 100048686

Hotplate		
Element(s)	Na	
SOP#		
Cocktail	40 mL 1:1 HCl	
sample(s)	weight (g)	volume (mL)
200128381_1	0.5031	50
200128381_2	0.5038	50
200128382_1	0.5024	50
200128382_2	0.5031	50

Date: 2026-06-22

Request ID: 100048686

LOI temp: 900 °C

sample	crucible (g)	crucible + sample (g)	after 900 °C	LOI (%)	correction
200128381	18.8059	20.8209	20.4176	20.01	0.7999

Na ICP analysis

sample	Dilution	conc. [mg/L]	Na_ppm [dried]	%Na [dried]
200128381_1	50/4	0.939	1457.64	0.15
200128381_2	50/4	0.927	1437.00	0.14

na result: 1447.32 ppm

na result: 0.14 %

na lot: *Manufacturer: Assurance, 27-159NAY exp: 2026-08-30*

Note(s):

Date: 2026-06-22

Request ID: 100048686

LOI temp: 900 °C

sample	crucible (g)	crucible + sample (g)	after 900 °C	LOI (%)	correction
200128382	16.8918	18.8976	18.4798	20.83	0.7917

Na ICP analysis

sample	Dilution	conc. [mg/L]	Na_ppm [dried]	%Na [dried]
200128382_1	50/3	0.906	1898.17	0.19
200128382_2	50/3	0.925	1934.23	0.19

na result: 1916.20 ppm

na result: 0.19 %

na lot: *Manufacturer: Assurance, 27-159NAY exp: 2026-08-30*

Note(s):

$$\text{ppm } M+ = [\text{conc.}][\text{volume}][\text{dilution}]/[\text{weight}]$$

$$\%M+ = [\text{conc.}][\text{volume}][\text{dilution}]/([\text{weight}] * 10,000)$$

$$\%MO = [\text{oxide factor}] * \%M+$$

A = crucible

B = crucible + sample

C = crucible + sample after drying

$$\%LOI = ([B-A] - [C-A]) / (B-A) * 100\%$$

$$\text{ppm } M+ \text{ [dried]} = \text{ppm } M+ / ([1 - (\%LOI)] / 100)$$

$$\%M+ \text{ [dried]} = \%M+ / ([1 - (\%LOI)] / 100)$$

$$\%MO \text{ [dried]} = \%MO / ([1 - (\%LOI)] / 100)$$