

Date: 2026-06-30

Request ID: 100048723

| Hotplate    |            |             |
|-------------|------------|-------------|
| Analyte(s)  | Cr6        |             |
| SOP#        |            |             |
| Cocktail    |            |             |
| sample(s)   | weight (g) | volume (mL) |
| 200128455_1 | 1.0025     |             |
| 200128455_2 | 1.0070     |             |
| 200128456_1 | 1.0031     |             |
| 200128456_2 | 1.0033     |             |

| Hotplate    |            |             |
|-------------|------------|-------------|
| Analyte(s)  | Cr2O3      |             |
| SOP#        |            |             |
| Cocktail    |            |             |
| sample(s)   | weight (g) | volume (mL) |
| 200128455_1 | 0.1023     |             |
| 200128455_2 | 0.1071     |             |
| 200128456_1 | 0.1043     |             |
| 200128456_2 | 0.1046     |             |

Date: 2026-06-30

Request ID: 100048723

Cr2O3 titration analysis

| sample      | weight (g) | titrant_volume (mL) | %Cr2O3 |
|-------------|------------|---------------------|--------|
| 200128455_1 | 0.1023     | 9.85297             | 22.278 |
| 200128455_2 | 0.1071     | 10.38171            | 22.421 |

FAS: 0.09125

cr2o3 result: 223495.00 ppm

cr2o3 result: 22.35 %

Cr6 titration analysis

| sample      | weight (g) | titrant_volume (mL) | %Cr6  |
|-------------|------------|---------------------|-------|
| 200128455_1 | 1.0025     | 5.37175             | 0.847 |
| 200128455_2 | 1.0070     | 5.39343             | 0.847 |

FAS: 0.09125

cr6 result: 8470.00 ppm

cr6 result: 0.85 %

Note(s):

Date: 2026-06-30

Request ID: 100048723

Cr2O3 titration analysis

| sample      | weight (g) | titrant_volume (mL) | %Cr2O3 |
|-------------|------------|---------------------|--------|
| 200128456_1 | 0.1043     | 10.35393            | 22.962 |
| 200128456_2 | 0.1046     | 10.39950            | 22.997 |

FAS: 0.09125

cr2o3 result: 229795.00 ppm

cr2o3 result: 22.98 %

Cr6 titration analysis

| sample      | weight (g) | titrant_volume (mL) | %Cr6  |
|-------------|------------|---------------------|-------|
| 200128456_1 | 1.0031     | 6.09530             | 0.961 |
| 200128456_2 | 1.0033     | 6.23328             | 0.982 |

FAS: 0.09125

cr6 result: 9715.00 ppm

cr6 result: 0.97 %

Note(s):

$$\%Cr(VI) = (1.733[mL\ FAS][N\ FAS])/([weight])$$

$$\%Cr_2O_3 = 1.462*[\%Cr(VI)]$$

$$\%Cr(III) = total\_Cr - Cr(VI)$$