

Date: 2026-07-14

Request ID: 100048713

Katanax		
Analyte(s)	Ni	
SOP#		
Acid Cocktail	30 mL 1:1 HCl	
sample(s)	weight (g)	volume (mL)
200128432_1	0.2585	250
200128432_2	0.2565	250
200128433_1	0.2574	250
200128433_2	0.2516	250
200128434_1	0.2599	250
200128434_2	0.2540	250

Date: 2026-07-14
Request ID: 100048713

Ni ICP analysis

sample	Dilution	conc. [mg/L]	Ni_ppm	%Ni
200128432_1	100/10	2.001	19352.03	1.94
200128432_2	100/10	2.084	20311.89	2.03

ni result: 19831.96 ppm

ni result: 1.98 %

ni lot: *Manufacturer: Assurance, 28-130NIT exp: 2026-10-30*

Note(s): *Retest: Nick reported a Ni analysis result of 2.07%.*

Date: 2026-07-14
Request ID: 100048713

Ni ICP analysis

sample	Dilution	conc. [mg/L]	Ni_ppm	%Ni
200128433_1	100/10	2.091	20308.86	2.03
200128433_2	100/10	2.027	20141.10	2.01

ni result: 20224.98 ppm

ni result: 2.02 %

ni lot: *Manufacturer: Assurance, 28-130NIT exp: 2026-10-30*

Note(s): *Retest: Nick reported a Ni analysis result of 2.10%.*

Date: 2026-07-14
Request ID: 100048713

Ni ICP analysis

sample	Dilution	conc. [mg/L]	Ni_ppm	%Ni
200128434_1	100/10	1.986	19103.50	1.91
200128434_2	100/10	1.935	19045.28	1.90

ni result: 19074.39 ppm

ni result: 1.91 %

ni lot: *Manufacturer: Assurance, 28-130NIT exp: 2026-10-30*

Note(s): *Retest: Nick reported a Ni analysis result of 1.85%.*

$$\text{ppm } M+ = [\text{conc.}][\text{volume}][\text{dilution}]/[\text{weight}]$$

$$\%M+ = [\text{conc.}][\text{volume}][\text{dilution}]/([\text{weight}]*10,000)$$

$$\%MO = [\text{oxide factor}]*\%M+$$