

Date: 2026-07-08
Request ID: 100048749

Katanax		
Analyte(s)	Ti	
SOP#		
Acid Cocktail	40 mL 1:1 HCl	
sample(s)	weight (g)	volume (mL)
200128532_1	0.2079	250
200128532_2	0.2087	250
200128533_1	0.2049	250
200128533_2	0.2049	250

Date: 2026-07-08
Request ID: 100048749

LOI temp: 900 °C

sample	crucible (g)	crucible + sample (g)	after 900 °C	LOI (%)	correction
200128532	16.8871	18.8943	18.5451	17.40	0.8260

Ti ICP analysis

sample	Dilution	conc. [mg/L]	Ti_ppm [dried]	%Ti [dried]
200128532_1	50/3	1.743	42290.01	4.23
200128532_2	50/3	1.848	44671.77	4.47

ti result: 43480.89 ppm

ti result: 4.35 %

ti lot: *Manufacturer: Assurance, 27-189TIT exp: 2026-08-30*

Note(s):

Date: 2026-07-08
Request ID: 100048749

LOI temp: 900 °C

sample	crucible (g)	crucible + sample (g)	after 900 °C	LOI (%)	correction
200128533	18.7362	20.7392	20.3627	18.80	0.8120

Ti ICP analysis

sample	Dilution	conc. [mg/L]	Ti_ppm [dried]	%Ti [dried]
200128533_1	50/3	1.873	46910.43	4.69
200128533_2	50/3	1.871	46841.56	4.68

ti result: 46876.00 ppm

ti result: 4.69 %

ti lot: *Manufacturer: Assurance, 27-189TIT exp: 2026-08-30*

Note(s):

$$\text{ppm } M+ = [\text{conc.}][\text{volume}][\text{dilution}]/[\text{weight}]$$

$$\%M+ = [\text{conc.}][\text{volume}][\text{dilution}]/([\text{weight}]*10,000)$$

$$\%MO = [\text{oxide factor}]*\%M+$$

A = crucible

B = crucible + sample

C = crucible + sample after drying

$$\%LOI = ([B-A]-[C-A])/(B-A)*100\%$$

$$\text{ppm } M+ \text{ [dried]} = \text{ppm } M+ / ([1-(\%LOI)/100])$$

$$\%M+ \text{ [dried]} = \%M+ / ([1-(\%LOI)/100])$$

$$\%MO \text{ [dried]} = \%MO / ([1-(\%LOI)/100])$$