

Date: 2026-08-03
Request ID: 100048851

Microwave				
Analyte(s)	Ni, Re			
SOP#				
Acid Cocktail	6 mL HNO3, 6 mL H2SO4			
Rack	7			
Vessel	Quartz			
Stir	90%			
Time (min)	Power (W)	T1 (°C)	T2 (°C)	P (bar)
15	1200	180	65	110
15	1500	240	65	110
30	1500	240	65	110
sample(s)	weight (g)	volume (mL)		
200128780_1	0.1030	200		
200128780_2	0.1010	200		

Ni, Re on carbon.

Date: 2026-08-03
Request ID: 100048851

Ni ICP analysis

sample	Dilution	conc. [mg/L]	Ni_ppm	%Ni
200128780_1	10/1	2.322	45087.38	4.51
200128780_2	10/1	2.278	45108.91	4.51

ni result: 45098.14 ppm

ni result: 4.51 %

ni lot: *Manufacturer: Assurance, 28-130NIT exp: 2026-10-30*

Re ICP analysis

sample	Dilution	conc. [mg/L]	Re_ppm	%Re
200128780_1	10/5	2.373	9214.24	0.92
200128780_2	10/5	2.269	8984.82	0.90

re result: 9099.53 ppm

re result: 0.91 %

re lot: *Manufacturer: Ricca, 840300A exp: 2027-04-30*

Note(s): *Ni, Re on carbon.*

$$\text{ppm } M+ = [\text{conc.}][\text{volume}][\text{dilution}]/[\text{weight}]$$

$$\%M+ = [\text{conc.}][\text{volume}][\text{dilution}]/([\text{weight}]*10,000)$$

$$\%MO = [\text{oxide factor}]*\%M+$$

Ni, Re on carbon.