

Date: 2026-08-06

Request ID: 100048866

Katanax		
Analyte(s)	Ni_total	
SOP#		
Acid Cocktail	30 mL conc. HCl	
sample(s)	weight (g)	volume (mL)
200128815_1	0.2066	250
200128815_2	0.2044	250

Hotplate		
Analyte(s)	Ni	
SOP#		
Cocktail	60 mL 1:1 HCl	
sample(s)	weight (g)	volume (mL)
200128815_1	0.2074	250
200128815_2	0.2037	250

Date: 2026-08-06
Request ID: 100048866

Ni ICP analysis

sample	Dilution	conc. [mg/L]	Ni_ppm	%Ni
200128815_1	50/1	1.703	102609.69	10.26
200128815_2	50/1	1.680	103072.33	10.31

ni result: #####

ni result: 10.28 %

ni lot: *Manufacturer: Assurance, 28-130NIT exp: 2026-10-30*

Ni_total ICP analysis

sample	Dilution	conc. [mg/L]	Ni_total_ppm	%Ni_total
200128815_1	50/1	1.730	104670.86	10.47
200128815_2	50/1	1.639	100232.39	10.02

ni_total result: #####

ni_total result: 10.25 %

ni_total lot: *Manufacturer: Assurance, 28-130NIT exp: 2026-10-30*

Note(s):

$$\text{ppm } M+ = [\text{conc.}][\text{volume}][\text{dilution}]/[\text{weight}]$$

$$\%M+ = [\text{conc.}][\text{volume}][\text{dilution}]/([\text{weight}]*10,000)$$

$$\%MO = [\text{oxide factor}]*\%M+$$